



The World's Most Compact Stair Lift

Get ahead of the curve by installing a Harmar Helix in your home today. Engineered to fit any curved staircase, the Harmar Helix utilizes True-Curve rail technology to precisely handle any twists and turns with a smooth, quiet operation. With its narrow profile and ability to be closely mounted to the wall, the Harmar Helix saves space and provides less intrusion to your family and guests..



RETRACTABLE SEAT BELT



NON-SKID FOOTREST



NARROW PROFILE



- **Narrow Profile:** Mounts 4.7" from wall. Low-profile floor mounts and a narrow chassis make this unit a space saving option, providing less intrusion on your stairway.
- **Built to Order:** Each Harmar Helix is custom-built to each location, ensuring a perfect fit every time.
- **Quicker Install:** Fewer floor bolts, simple connecting rails, and a light-weight chassis allow for easy installation and minimal disruptions for you.
- **Color Options:** Every lift comes with the ability to choose custom colors and fabric options.
- **Made in the USA.**

SPECIFICATIONS	HELIX
Weight Capacity	350 lbs.
Travel	16', 18', and 24' Standard (depending on configuration); Extra rail up to 164'
Control in Armrest	Yes (left or right hand operation)
Folded Width	16" from wall
Footrest Size	11" x 12"
Armrests	Fully flip back, folded armrests
Floor to Seat Height	25"
Seat Size	9" W x 15" D x 20.5" H
Operation Power	24V DC Battery
Speed	20 fpm
Incline Limits	Up to 60°
Outdoor	Yes
Drive System	Steel Rack and Pinion
Safety Design	Complies with ASME A18.1, CAN/CSA B44.1, ASME A17.5
Safety Features	<ul style="list-style-type: none"> • Direction Limit Switches • Final Limit Switch • Footrest Obstruction Sensors • Seat Swivel Interlock • Constant Pressure Control • Retractable Seat Belt • Chassis Obstruction Sensors • Overspeed Brake
Warranty	<ul style="list-style-type: none"> • Limited Lifetime: Drive System, Structure • 2-Years: Parts • 1-Year: Batteries



1 Year Limited Warranty
on New Equipment Installs

An independently owned and operated franchise business.

- Stair Lifts
- Vertical Platform Lifts
- Pool Lifts
- Patient Lifts
- Ramps
- Auto Lifts
- Bath Modifications
- Commercial
- Residential
- And More